



**ken-a-vision**<sup>®</sup>  
KNOWLEDGE THROUGH VISION

## Megapixel Vision Viewer 7890UM

The Megapixel Vision Viewer 7890UM has a high resolution camera and a USB 2.0 connection for a fast connection to your computer or laptop. It comes complete with Applied Vision™ 4 Software. AV4is MSI compatible and offers image capture, movies, time-lapse, rotate, brightness, contrast and much more.



### Features and Benefits

- **Rubber Knurled Focusing System** Better grip for smoother focusing
- **Pre-Attached Cabling** No cable to lose
- **Resolution** 1280 x 1024 SXGA
- **Ergonomically Designed Swivel Head** For ease of positioning
- **Lens** 6mm glass
- **Applied Vision™ Software** Full Screen viewing

#### 7890UM SPECIFICATIONS

Weight	4 lbs / 1.82 kgs	Connection	Pre-attached 5' (1.52 m) USB Cable
Base Dimension	7" x 7" (18 cm x 18 cm)		USB 2.0 / 1.1
Lens	6mm Glass	Exposed Cable Length	3' 7" (1.12m) USB Cable
Neck Length	20" (51 cm)	Real-Time Video	Up to 30 Frames per second
Focal distance	8 mm to infinity	Microscope Adapter	28mm included
Head Type	Swivel	Software	Applied Vision™
Resolution	1280 x 1024 SXGA	Warranty	3 Years

Specifications subject to change without notice



**ken-a-vision**<sup>®</sup>  
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: [info@ken-a-vision.com](mailto:info@ken-a-vision.com)