

# FlexCam by Ken-A-Vision®

## FlexCam iCam Digital 910-171-102



The FlexCam iCam Digital features both S-Video and USB outputs for easy connection to a PC, Macintosh, television, VCR or LCD projector. This innovative, low-cost camera is perfect for presenting objects, images, text, sending video email and creating video portfolios. It includes powerful imaging software to facilitate video and still image capture and editing. The FlexCam iCam is ideal for viewing microscope images. The integrated focus knob is sized and shaped perfectly for use with the included microscope adapters.

### Additional Features and Benefits

- PC and Mac Compatible
- USB Digital-output
- Ease of Connection to Most Microscopes
- Superior Positioning

#### 910-171-102 SPECIFICATIONS

Weight	4 lbs / 1.81 kg	Resolution	330 Lines
Dimension	7" x 7.5" x 20"	Real-Time Video	up to 30 Frames per second
	18 cm x 19 cm x 50.8 cm	Output Signal	S-Video
Gooseneck Length	20" (50.8 cm)		Digital (USB2.0)
Lens	8 mm D-mount	Software	ArcSoft VideoImpression
Illumination	3 Lux		PhotoStudio
Image Type	CMOS	Microscope Adapters	26 mm, 28 mm & 34 mm
White Balance	Auto	Warranty	2 Years
Grain Control	Auto		
Iris	Auto with flicker-free performance		

Specifications subject to change without notice



**ken-a-vision®**  
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: [info@ken-a-vision.com](mailto:info@ken-a-vision.com)

kav.fly.910171102.107v3