



ken-a-vision®
KNOWLEDGE THROUGH VISION

Cordless PrepScope

T-1201C

This compact cordless microscope is designed with optical quality and ruggedness to withstand student use. The Cordless PrepScope offers ease of use and provides solutions for classrooms with a lack of electrical outlets. It is also great for use in your outdoor lab.



Features and Benefits

- **Prism Head**
 - **Head**
 - **Rechargeable**
 - **Slip Clutch**
 - **High density plastic covered base**
 - **Power**
- For true optical clarity (no mirror)
Monocular
NiMH batteries
Prevents damage to focusing mechanism
Protects the metal base from chipping
40 hours for a single charge or use with power supply for unlimited time



T-1201C SPECIFICATIONS

Weight	4 lbs / 2 Kg	Light Source	Ken-A-Vision Cool Lighting
Size	12"Hx 6"D x 4"W	Light Board	100,000 hours
Prism Head	45° angle, rotates 360°	Power Source	Rechargeable NiMH batteries
Head	Monocular	Power Source Life	500 cycles (approximately 3 years)
Eyepiece	10x Widefield	Base	High density plastic covered metal
Objective Lens	4x, 10x, 40xR DIN	Warranty	10 Years
Stage Size	90 mm x 90 mm		
Coaxial Focusing	One Knob-coarse & fine		
Diaphragm	6 aperture disc		
Stage	Spring loaded stage clips		

Specifications subject to change without notice



ken-a-vision®
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com