



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

Cordless Dual Purpose Comprehensive Scope T-1933C

The new Cordless Comprehensive Dual Purpose Scope combines a standard compound microscope with limited stereo microscope capability for one low price.



Features and Benefits

- Separate Top Light for Viewing Large, Macro or Opaque Specimens
- Bottom Light for Standard Microscope Slides
- Light Rheostat for Better Light Control
- Seidentöpf Binocular Head



T-1933C SPECIFICATIONS

Weight	15 lbs / 8 Kg	Stage	x/y Floating
Base Material	All metal construction	Above Stage Lighting	Yes
Dimension	8.5"L X 6.5"W X 15"H	Below Stage Lighting	Yes
	21.6 cm X 16.5 cm X 38.1 cm	Light Source	Ken-A-Vision Cool Lighting
Head Type	Seidentöpf Binocular Head	Bulb Life	100,000 hours
Eyepiece	10x Widefield	Power Source	Rechargeable NiMH Batteries
Objective Lens	2x, 4x, 10x, 40xR DIN	Warranty	10 Years, Excluding Batteries & Lamp
Focal Adjustment	Coaxial Coarse & Fine		

Specifications subject to change without notice



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com