



ken-a-vision®
KNOWLEDGE THROUGH VISION

The Professor

ESH101 & ESH111

The Professor by Ken-A-Vision is the world's first microscope designed with above light for viewing opaque specimens. The Professor is a low cost solution for introducing students to the world of microscopy. With a compact size, battery power and glass optics, The Professor by Ken-A-Vision is the ideal microscope for working in the classroom or in the field. The Professor is available for purchase as a single unit or in packs of ten.



Additional Features and Benefits

- Top Light for Micro-magnifications
- High Impact Plastic Construction
- 4x, 10x, 40x Glass objective lens
- 10x Glass eyepiece
- 300 Hours of operation utilizing three AA batteries (included)
- Dust cover
- 10 Safety Slides (5 prepared, 5 blank)

PROFESSOR ESH101 SPECIFICATIONS

Weight	1.25 lbs / 0.57 kg	Above Stage Lighting	Ken-A-Vision Cool Lighting
Dimension	11.75"H X 4.75"W X 6"	Illuminator - Sub Stage	Ken-A-Vision Cool Lighting
	30 cm X 12 cm X 15.2 cm	Power Source	3 AA Batteries (included)
Eyepiece	10x	ESH101	Single Unit
Optics	Glass	ESH111	10 Microscopes
Head	45° angle single fixed	Accessories	10 Safety Slides (5 prepared, 5 blank)
Objective Lens	4x, 10x, 40x	Warranty	1 Year
Stage Size	3.25" X 3.5" / 8.5 cm X 9.0 cm		
Diaphragm	6 aperture disc		

Specifications subject to change without notice



ken-a-vision®
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com

kav.fly.ESH101.107v4