



ken-a-vision
KNOWLEDGE THROUGH VISION

Video Flex[®]

2100 Explorer

The 2100 Explorer will meet all your needs without exceeding your budget. It's the ideal solution when seeking portability and versatility. The 40-inch neck can be coiled into a base for microscope applications, or can be wrapped around a telescope, object, or person to provide a unique point of view.



Additional Features and Benefits

- **Baseless Construction** Unlimited positioning and wrapping around objects
- **Light Weight** Transports easily for remote video capture
- **Poly-Carbonate Build** Allows for rugged conditions and environment
- **C-Mount Lens** Interchangeable lens system accommodates all c-mount lenses
- **Quick Focus™** 1/4 turn moves from 6mm to infinity
- **450 Lines of Resolution** For crisp clear pictures
- **Portable** At 1.5 pounds travel and explore

2100 SPECIFICATIONS

Weight	1.54 lbs / 0.7kg	Resolution	450 Lines
Base	Neck is the base	Magnification	50:1 and Greater
Neck Length	40 " / 102 cm	Exposure	Fixed Contrast
Head Type	90 degree Fixed	Output Signal	Composite Video
Lens	8 mm	Eyepiece Adapters	28mm & Built-in 34.5 mm
Focus	¼ turn Quick Focus™	Cable	12' Video
Focal distance	¼" / 6mm to Infinity	Warranty	2 Years
Light Sensitivity	1.5 Lux		

Specifications subject to change without notice



ken-a-vision
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com